

Introducing the next generation in Airgel. Our iQ® Blue Nasal Vented and Non-Vented masks feature our patented AIR°gel with *Advanced Cushion Technology* and flexible shell with *Custom Fit Technology*™. Together, they provide an optimal fit with maximum comfort. Unlike masks with rigid shells, iQ Blue adjusts to truly fit the contours of your patients' face.

The Smart Choice for a Custom Fit

Custom Fit Technology

The shapeable shell of each mask can be easily molded to fit individual needs.
Flex, fit, forget.



Advanced Cushion Technology (A.C.T.)

Our ergonomic cushion technology delivers an effective seal for the most comfortable fit.

360° Rotating Elbow

For optimal tube positioning and easier mobility.

Clearest Field of Vision Possible

No forehead pad required.

AIR°gel

Maximum comfort, minimum weight. Our patented AIR°gel minimizes pressure points and skin irritation.

Quick-Release Headgear

Lightweight, breathable headgear slides off easily with one hand.

Vented and Non-Vented Versions Available



What makes our iQ Blue Mask different?

- *Custom Fit Technology* is unique to Sleepnet®. Our soft-shell mask has a shapeable wire allowing for optimal flexibility, fit, and comfort.
- Our *AIR°gel with Advanced Cushion Technology* creates an ergonomic gel cushion that reduces pressure points on the skin.
- Clearest field of vision possible — users can even wear glasses with this mask.
- No forehead pad, yet retains outstanding stability.

About Advanced Cushion Technology (A.C.T.)

Our patented AIR°gel with an ergonomic streamlined design and luxurious feel reduces pressure while maintaining an effective seal.



How does iQ Blue stand apart from the competition?

| Features | Sleepnet iQ Mask | ResMed Mirage™ Softgel Mask | Fisher & Paykel Forma Mask | Philips Respironics Comfort Gel Mask |
|-----------------------------|------------------|-----------------------------|----------------------------|--------------------------------------|
| Custom Fit Technology | ● | | | |
| Advanced Cushion Technology | ● | | | |
| Minimal forehead contact | ● | | | |
| Easy universal sizing | ● | ● | ● | ● |
| Breathable headgear | ● | | | |
| Patented AIR°gel | ● | | | |
| 100% made in the USA | ● | | | |

iQ Blue Mask Specifications:

| Product Number | Mask Style | Headgear Style | Mask Size | Tube Length | Mask Weight | Mask Dimensions |
|----------------|------------|----------------|-----------|-------------|-------------|--------------------------|
| 50647 | Vented | 3-Point | OS | 13" | 63.8 g | 3.04"D x 2.95"W x 2.60"H |
| 50655 | Vented | StableFit | OS | 13" | 63.8 g | 3.04"D x 2.95"W x 2.60"H |
| 50250 | Vented | 3-Point | OS | 6" | 63.8 g | 3.04"D x 2.95"W x 2.60"H |
| 50255 | Vented | StableFit | OS | 6" | 63.8 g | 3.04"D x 2.95"W x 2.60"H |
| 50660 | Non-Vented | 3-Point | OS | 6" | 63.9 g | 3.08"D x 3.13"W x 2.50"H |



Mojo®

Your Mask, Your Way

Our full-face mask also features AIR°gel.



Contact us today.

To place your order, call an authorized distributor.

Call 1.800.742.3646

Visit www.sleepnetmasks.com to learn more.

Sleepnet Corporation, 5 Merrill Industrial Drive, Hampton, New Hampshire 03842 USA



Scan code to view iQ Fitting Tips video.